



ShanghAI Lecture Series

15 November 2012, 10:00 AM CET

# Design Principles for Intelligent Rehabilitation Robots

**Robert Riener**

Sensory-Motor Systems Lab

Institute of Robotics and Intelligent Systems, ETH Zurich  
University Hospital Balgrist, University of Zurich

# Disadvantage of Manual Training



For the Therapist

- Physically exhausting
- Ergonomically inconvenient



For the Patient

- Limited training duration
- Gait pattern not optimal



SCI Center, Balgrist

# Human-Robot Gait Rehabilitation




# Gait Rehabilitation "Robots"



GaitTrainer



G-EO



Haptic Walker



HAL




Auto  
ambulator

# ARMin III



## Exoskeletal Robot with 7 Degrees of Freedom



ETH Zurich/Balgrist, Hocoma AG

Nef, Riener et al. 2006-2011

# Arm Rehabilitation Robots



MIT Manus



Bi-Manu-Track



MGA



Salford PMA



Haptic Master, GENTLE/s



PERCRO Exoskeleton




KIST Arm



# Sensory-Motor Interaction



Human



Machine

Interaction

**Focus: Exoskeletons for Clinical Rehabilitation**  
**Setups: Lokomat, ARMin**

# Sensory-Motor Characteristics



## To Be Known about Human Physiology




- Nervous system is plastic

# Sensory-Motor Characteristics



## To Be Known about Human Physiology




- Nervous system is plastic

# Sensory-Motor Characteristics



## To Be Known about Human Physiology



- Nervous system is plastic
- Participation increases plasticity

# Passive vs Active Training



## Limitations of Passive Training of Healthy Subjects

- Physical guidance hinders motor learning of walking balance

Domingo &

Ferris, 2009


- High frequency guiding is detrimental for learning of arm movements (guidance hypothesis)

Winstein et al. 1994; Marchal-Crespo & Reinkensmeyer 2008


## Assist-as-Needed Training (AAN)

- AAN shows higher level of recovery step number, periodicity, and consistency (27 mice)

Cai et al. 2006



# Lokomat: Patient-Cooperative Control



## Path Control


- Robot behaves **assistive**, **corrective** or **transparent**, when needed
- Free timing for patient
- Support patient, but do not restrict patient

Duschau-Wicke, Vallery, Riener 2009, 2010, 2011

# Path Control Increases Participation



## Muscle Activity



11 incomplete SCI subjects

## Heart Rate

Relative increase of heart rate


Position  
control

Path  
control

# Sensory-Motor Characteristics



## To Be Known about Human Physiology




- Nervous system is plastic
- Participation increases plasticity
- Repetition without repetition

# Path Control Enhances Variability




**“Repetition without Repetition”**

Position Control




Path Control



# Sensory-Motor Characteristics



## To Be Known about Human Physiology




- Nervous system is plastic
- Participation increases plasticity
- Repetition without repetition
- Provide large ROM, avoid joint stress

# Main Problem with Exoskeletons



## Alignment of Robotic and Human Joints




Human

Robot


# Main Problem with Exoskeletons



## Alignment of Robotic and Human Joints



Human  
Robot




Human  
Robot


# Main Problem with Exoskeletons



## Alignment of Robotic and Human Joints




Human  
Robot




Human  
Robot

# Humerus Motion




# Humerus Motion





Body height: 1.7 m


Center of the Glenohumeral Joint (Caput Humeri)




# ARMin II Shoulder Kinematics




# ARMin III: Novel Shoulder Actuation



# ARMin III: Novel Shoulder Actuation




# ARMin III: Novel Shoulder Actuation



# Sensory-Motor Characteristics





## To Be Known about Human Physiology



- Nervous system is plastic
- Participation increases plasticity
- Repetition without repetition
- Provide large ROM, avoid joint stress
- Each patient is different




# ARMin III: Adjustments



# Lokomat: Adaptive Force Field




Adapt Field with a Iterative Learning Controller



# Sensory-Motor Characteristics




## To Be Known about Human Physiology



- Nervous system is plastic
- Participation increases plasticity
- Repetition without repetition
- Provide large ROM, avoid joint stress
- Each patient is different
- Motivation increases plasticity

# Human-Robot Cooperation



# Sensory-Motor Characteristics



## To Be Known about Human Physiology




- Nervous system is plastic
- Participation increases plasticity
- Repetition without repetition
- Provide large ROM, avoid joint stress
- Each patient is different
- Motivation increases plasticity
- Train activities of daily living



# Task Specific Training



**Chronic Stroke Patient (FMA=26)**



# Task Specific Training



# Conclusion and Outlook



## Taking into Account Human Physiology

- Robots can allow efficient & intensive & individual training
- Robot can cooperate to keep the patient active
- Robots can motivate the patient

## Chances

- Improve health status and quality of life (patients)
- Reduce workload of clinical staff (therapists, nurses)

## To do

- Clinical evaluation studies on patients